


Enhanced Development and Control of Continuous Processes



METTLER TOLEDO



Agenda

Continuous Flow Chemistry - Analysis Challenges

- FlowIR™ for Continuous Flow Chemistry
- Application One - Oxazole Formation
- Application Two - Fluorination Reaction
- Application Three - Accurate Control of Reagent Addition in Multi-step Segmented Flow Processes Using Inline Infra-red Monitoring
- Summary

Continuous Chemistry - Analysis Challenges

One of the main analytical challenges is real-time inline monitoring for enhanced understanding and control of reactions, especially when dealing with multistep sequences



- Chemical
 - Information concerning **reaction monitoring** is highly desirable, as this normally requires taking a sample offline for TLC, LCMS, UV, etc.
 - Inline monitoring would greatly assist with rapid optimization procedures
- Chemical
 - Information concerning the **stability of reactive intermediates** is highly desirable, as these cannot always be detected by other methods
- Technical
 - Information concerning **dispersion and diffusion** is highly desirable as these are unavoidable side effects of continuous flow - very important to characterize these parameters for multistep reaction sequences

Agenda

- Continuous Flow Chemistry - Analysis Challenges

FlowIR™ for Continuous Flow Chemistry

- Application One - Oxazole Formation
- Application Two - Fluorination Reaction
- Application Three - Accurate Control of Reagent Addition in Multi-step Segmented Flow Processes Using Inline Infra-red Monitoring
- Summary

FlowIR™ for Flow Chemistry



FlowIR™

- Designed specifically for flow chemistry
- High performance with no utility requirement (except power)
- Universal (fits with any flow system)

- Determination of dispersion effects
 - Used for stoichiometric control
- Multi-step/stage reaction control
- Faster structural information in real time
- Immediate determination of steady-state conditions
- Immediate detection of reaction start/end point
- Information on reaction mechanism/pathway
- Detection of transient intermediates
- Information without the need to sample



Agenda


- Continuous Flow Chemistry - Analysis Challenges
- FlowIR™ for Continuous Flow Chemistry

Application One - Oxazole Formation


- Application Two - Fluorination Reaction
- Application Three - Accurate Control of Reagent Addition in Multi-step Segmented Flow Processes Using Inline Infra-red Monitoring
- Summary

Oxazole Formation

Monitoring a continuous oxazole formation within a microfluidic device (10 μ L)



- Oxazole formation is monitored by following its unique carbonyl absorbance at 1721cm⁻¹
- Benzylamine column shown to be an effective scavenger of excess acid chloride (used to drive reaction to completion)
- Dispersion effect caused by columns clearly shown in product trend curve



Agenda


- Continuous Flow Chemistry - Analysis Challenges
- FlowIR™ for Continuous Flow Chemistry
- Application One - Oxazole Formation

Application Two - Fluorination Reaction

- Application Three - Accurate Control of Reagent Addition in Multi-step Segmented Flow Processes Using Inline Infra-red Monitoring
- Summary

Fluorination Reaction Using DAST

Experimental setup



CaCO₃ and SiO₂ column provides an inline quench and cleanup

ReactIR™ 45m



Vapourtec R2+/R4 system



Convection Flow Coil (CFC) reactor

DS Micro Flow cell

Reaction output

Fluorination Reaction

Dispersion effects are clearly monitored



- Expected time for product to come out of the reaction line: 20 mins (A), [4mL total volume at 0.2mL/min = 20 mins]
- Actual 36 mins (A+B)
- Additional time B is required due to dispersion, diffusion and chromatographic effects caused by polymer-supported reagents

Agenda


- Continuous Flow Chemistry - Analysis Challenges
- FlowIR™ for Continuous Flow Chemistry
- Application One - Oxazole Formation
- Application Two - Fluorination Reaction

Application Three - Accurate Control of Reagent Addition in Multi-step Segmented Flow Processes Using Inline Infrared Monitoring

- Summary

Accurate Control of Reagent Addition in a Multi-step Process


- The unavoidable dispersion of the reaction “plug” is a significant issue when performing multi-step sequences in flow
- The controlled addition of exact stoichiometries of reagents to a product stream is therefore challenging with the current commercially available flow equipment
- Poor control is wasteful on small scale (chiral/expense/toxic material used in excess) and requires additional purification.



- Current situation - use the IR to generate a dispersion curve of the intermediate and manually switch the pump on

Accurate Control of Reagent Addition in a Multi-step Process

- Alternative: use real-time concentration of a product specific band from the ReactIR™ and convert it into a flow rate that allows the third stream to be dispensed proportional to the concentration of the intermediate via LabVIEW



Can We Add a Third Stream with Accurate 1:1 Stoichiometry?

X = stoichiometry between the two components

F = flow rate of output process


[D] = conc of third stream component

A = peak height value measured by IR

k = conversion factor

$$\text{Flow rate of third stream} = \frac{XFk}{[D]} \cdot A$$

- K determined from a concentration screen of the intermediate
- Test system: two inert compounds (4-chlorobenzophenone & 3-methyl-4-nitroanisole)






$$\text{FR} = ((1 \times 0.25 \times 2210) / 0.0833) \times A$$
$$= 6630 \times A$$

Therefore labview will multiply the live peak height value of 4-chlorobenzophenone by 6630 to generate a concentration dependent flow rate in $\mu\text{L}/\text{min}$


Can We Add a Third Stream with Accurate 1:1 Stoichiometry?

4-chlorobenzophenone


3-methyl-4-nitroanisole



- Third stream successfully added with 1:1 stoichiometry for >97% of the material
- Limitation at the ends of the dispersion curves arises from inaccuracy of piston pumps at very low flow rates
- Absolute ratios determined using ¹H NMR spectroscopy to confirm results




Application to yne-one Project



No IR control


- 10 equiv toxic hydrazine used
- Visual observation used to manually switch the third pump
- Extensive purification required

With IR control



- Amount of toxic hydrazine pumped reduced to 3 equiv.
- Reduced reaction temperature to avoid polymerisation of terminal acetylene
- Plug of silica gel allows chromatographic separation with IR detection
- Crystalline material obtained from both processes

Application to yne-one project



Agenda

- Continuous Flow Chemistry - Analysis Challenges
- FlowIR™ for Continuous Flow Chemistry
- Application One - Oxazole Formation
- Application Two - Fluorination Reaction
- Application Three - Accurate Control of Reagent Addition in Multi-step Segmented Flow Processes Using Inline Infra-red Monitoring

Summary

Summary

- ReactIR™ DS Micro Flow Cell (FlowIR™) has been shown to instantaneously obtain highly molecular-specific chemical information. This can provide a unique insight into continuous flow chemistry
- The case studies demonstrate that FlowIR™ can be used with a range of different scale flow reactors:
 - Microscale - 10µL (Future Chemistry)
 - Meso scale flow reactors (Uniqsis, Vapourtec)
 - Large kilo lab flow reactors (Alfa Laval)
- The ability to track dispersion curves now makes it feasible to accurately apply a third stream via a feedback controlled flow that exactly matches the required amount of reagent to the output of the first flow reactor. This leads to improved product yield and reduced purification steps.