

Miniaturization of a Luciferase Reporter Gene Assay Shows Enhanced Assay Performance With Considerable Cost Savings



SnAPPShots

*A brief report
from the Corning
Applications Group*

David Hayes
Mark E. Rothenberg
Michael Briggs
Todd Upton


Corning Incorporated
Life Sciences

Introduction and Purpose

Assay miniaturization has become an important goal for scientists attempting to reduce costs associated with drug development. By moving to a Low Volume (LV) format one effectively reduces the waste associated with Normal Volume (NV) assays. Reducing assay volumes from 100 μL to 10 μL , like those found in 96 and low volume 384 well microplates respectively, can drop reagent costs linearly with volume (1). The purpose of this study is to demonstrate how miniaturizing a Luciferase reporter gene assay from a traditional 45 μL /well down to 17.5 μL /well utilizing the Corning® 384 well LV microplates impacts assay performance.

Methods and Results

Using the Corning 384 LV well tissue culture (TC)-treated flat bottom microplates (Cat No. 3826), HEK 293 CRE-luc cells (Panomics Cat. No. RC0007) were seeded at a density of 2,500 cells/well in either 10 (LV) or 20 μL (NV). The cells were cultured in IMDM without phenol red with 10% FBS and 0.1 mg/mL hygromycin. This cell line is stably transfected with a luciferin reporter gene driven by a promoter responsive to the cyclic AMP response element-binding (CREB) protein. The cells were induced with 20 μM forskolin in a volume of 2.5 or 5 μL by incubation overnight, in a humidity controlled incubator at 37°C with 5% CO_2 . Prior to reading, the plates were equilibrated to room temperature for 30 minutes, after which 5 or 20 μL of SteadyLite HTS Luciferase reagent (Perkin Elmer™) were added to the appropriate wells and the plates read on a Molecular Devices Aquest® plate reader.



Corning 384 Well Low Volume Flat Bottom Microplate (Cat. No. 3826)

By miniaturizing the assay we observed a considerable cost savings in terms of reagent and compound use (refer to Table 1). The cost savings are realized by the 75% reduction in the SteadyLite reagent utilized in the LV assay which equates to a cost/well that is 24% of a NV plate. There was also minimal impact on the Z' score (Figure 1) and similar, if not higher, signal to background ratio in the LV versus NV assay plates (Figure 2). The geometry of the wells of Corning's 384 well tissue culture (TC)-treated flat bottom microplates functions not only to allow cost savings but also to concentrate the signal in the LV assay and improve the S:B ratio.

CORNING

Figure 1. Z' Score





Figure 2. S:B



Conclusions

- ▶ Miniaturization leads to a 75% reduction in SteadyLite usage as well as a considerable cost savings using the low volume assay format.
- ▶ Assay robustness in LV plates is equal to or better than that seen in NV plates.
- ▶ The overall geometry and LV format equate to both cost savings and enhanced assay performance.

Reference

1. Tina K. Garyantes, 1536-well assay plates: when do they make sense? DDT Vol. 7, No. 9 May 2002 editorial.

Table 1.

	Normal Volume Assay	Low Volume Assay
mL of SteadyLite/kit	1000	1000
µL SteadyLite/well	20	5
Number of wells/1000 mL	50000	200000
List price of reagent kit	\$2,688	\$2,688
Cost per well	0.054	0.013

Corning has been transforming possibilities into breakthrough realities for over 150 years.

CORNING

Corning Incorporated Life Sciences

Tower 2, 4th Floor
900 Chelmsford St.
Lowell, MA 01851
t 800.492.1110
t 978.442.2200
f 978.442.2476

www.corning.com/lifesciences

Worldwide Support Offices

ASIA / PACIFIC

Australia
t 61 2-9416-0492
f 61 2-9416-0493

China
t 86 21-3222-4666
f 86 21-6288-1575

Hong Kong
t 852-2807-2723
f 852-2807-2152

India
t 91-124-235 7850
f 91-124-401 0207

Japan
t 81 (0) 3-3586
1996/1997
f 81 (0) 3-3586
1291/1292

Korea
t 82 2-796-9500
f 82 2-796-9300

Singapore
t 65 6733-6511
f 65 6861-2913

Taiwan
t 886 2-2716-0338
f 886 2-2716-0339

EUROPE

France
t 0800 916 882
f 0800 918 636

Germany
t 0800 101 1153
f 0800 101 2427

The Netherlands
t 31 20 655 79 28
f 31 20 659 76 73

United Kingdom
t 0800 376 8660
f 0800 279 1117

All Other European Countries
t 31 (0) 20 659 60 51
f 31 (0) 20 659 76 73

LATIN AMERICA

Brasil
t (55-11) 3089-7419
f (55-11) 3167-0700

Mexico
t (52-81) 8158-8400
f (52-81) 8313-8589