

CORNING® 96 WELL MICROPLATE SELECTION CHART

Page Number	AVAILABLE SURFACES
Not Treated	
High Bind	
Medium Bind	
NBS™	
Carbo-BIND™	
Sulphydryl-BIND™	
Universal-BIND™	
DNA-BIND®	
Amine	
TC-Treated	
Corning CellBIND® Surface	
Ultra-Low Attachment Surface	
Poly-D-Lysine	
UV Transparent	

96 Well Clear Polystyrene Microplates

[illegible]

96 Well Solid Black Polystyrene Microplates

Round Bottom	9	■												
Flat Bottom	9	■	■		■					■				
Flat Bottom Half Area	9	■			■					■				

96 Well Black Clear Bottom Polystyrene Microplates

Flat Bottom	9-10	■	■		■					■	■		■	
Flat Bottom Half Area	10	■			■					■			■	

96 Well Solid White Polystyrene Microplates

Round Bottom	9	■			■									
Flat Bottom	9	■	■		■					■				
Flat Bottom Half Area	9	■			■					■				

96 Well White Clear Bottom Polystyrene Microplates

Flat Bottom	9-10	■			■					■			■	
Flat Bottom Half Area	10	■			■					■				

96 Well Stripwell Polystyrene Microplates

[illegible]

96 Well UV Microplates

[illegible]

96 Well Clear Flexible Polyvinyl (PVC) Microplates

[illegible]

96 Well Polypropylene Microplates

[illegible]

96 Well Polypropylene Storage Blocks

[illegible]

CORNING® 384 AND 1536 WELL MICROPLATE SELECTION CHART

		AVAILABLE SURFACES								
	Page Number	Not Treated	High Bind	Medium Bind	NBS™	TC-Treated	Corning CellBIND® Surface	Ultra-Low Attachment Surface	Poly-D-Lysine	UV Transparent
384 Well Clear Polystyrene Microplates										
Flat Bottom	16	■	■		■	■			■	
384 Well Solid Black Polystyrene Microplates										
Flat Bottom	16	■			■	■				
Round Bottom, Low Volume	16	■	■	■	■					
Flat Bottom, Low Volume	16	■			■	■				
384 Well Black Clear Bottom Polystyrene Microplates										
Flat Bottom	17	■			■	■	■	■	■	
Flat Bottom, Low Volume	17	■			■	■				
384 Well Solid White Polystyrene Microplates										
Flat Bottom	16	■			■	■				
Round Bottom, Low Volume	16	■		■	■					
Flat Bottom, Low Volume	16				■	■				
384 Well White Clear Bottom Polystyrene Microplates										
Flat Bottom	17	■			■	■			■	
384 Well UV Microplates										
Flat Bottom	17									■
384 Well Polypropylene Microplate										
Round Bottom	18	■								
Conical Bottom	18	■								
384 Well Polypropylene Storage Block										
Round Bottom	18	■								
V-Bottom	18	■								
1536 Well Clear COC Echo Qualified Microplates										
Flat Bottom	20	■								
1536 Well Solid Black Polystyrene Microplates										
Flat Bottom	19	■			■	■				
Round Bottom	19	■								
1536 Well Black Clear Bottom Polystyrene Microplates										
Flat Bottom, High Base	20	■			■	■	■			
Flat Bottom, Low Base	20	■				■	■			
1536 Well Solid White Polystyrene Microplates										
Flat Bottom	19	■			■	■	■			
Round Bottom	19	■								