

Corning® Mini Bioreactor

CORNING

The new Corning Mini Bioreactors are ideal for high throughput process optimization for suspension cell culture. The product consists of the Corning 50 mL Centrifuge Tube with a vented cap. Use the product in the following applications:

- ▶ Cell line development
- ▶ Clone selection
- ▶ Media optimization
- ▶ Recombinant proteins

Features and Benefits


- ▶ Polypropylene tube with large marking spot to clearly identify tube contents and experimental parameters
- ▶ Polypropylene cap with 4 vents; hydrophobic membrane provides gas exchange
- ▶ Sterile, nonpyrogenic and RNase-/DNase-free

New!



Corning® Mini Bioreactor Ordering Information

Cat. No.	Description	Qty/Cs
431270	Mini Bioreactor, 50 mL with hydrophobic vented cap	300



www.corning.com/lifesciences/solutions

At Corning, cells are in our culture. In our continuous efforts to improve efficiencies and develop new tools and technologies for life science researchers, we have scientists working in Corning R&D labs across the globe, doing what you do every day. From seeding starter cultures to expanding cells for assays, our technical experts understand your challenges and your increased need for more reliable cells and cellular material.

It is this expertise, plus a 160-year history of Corning innovation and manufacturing excellence, that puts us in a unique position to offer a beginning-to-end portfolio of high-quality, reliable cell culture consumables.

For additional Corning product, technical, or distributor information, please e-mail us at cctech@corning.com, visit our website www.corning.com/lifesciences or call + 31 20 659 6051 or contact your local Corning sales office listed below.

Corning Incorporated Life Sciences

Corning BV
Fogostaat 12
1060 LJ Amsterdam
The Netherlands
Phone: +31 0 20 659 60 51
Fax: +31 20 659 76 73
cceurnl@corning.com
www.corning.com/lifesciences

Worldwide Support Offices

ASIA/PACIFIC

Australia/New Zealand
t 0402-794-347

China
t 86 21 2215 2888
f 86 21 6215 2988

India
t 91 124 4604000
f 91 124 4604099

Japan
t 81 3-3586 1996
f 81 3-3586 1291

Korea
t 82 2-796-9500
f 82 2-796-9300

Singapore
t 65 6733-6511
f 65 6861-2913

Taiwan
t 886 2-2716-0338
f 886 2-2516-7500

EUROPE

France
t 0800 916 882
f 0800 918 636

Germany
t 0800 101 1153
f 0800 101 2427

The Netherlands
t 31 20 655 79 28
f 31 20 659 76 73

United Kingdom
t 0800 376 8660
f 0800 279 1117

All Other European Countries

t 31 (0) 20 659 60 51
f 31 (0) 20 659 76 73

LATIN AMERICA

Brasil
t (55-11) 3089-7419
f (55-11) 3167-0700

Mexico
t (52-81) 8158-8400
f (52-81) 8313-8589

The Corning Family of Brands

CORNING

costar®

PYREX®

AXYGEN™

GOSSELIN™